

COMPARISON TABLE

COMPARED PRODUCTS	M-BASIC	M-DIAGNOSTIC	M-SCIENCE
			
POWER SUPPLY	230V 50/60Hz; 120V 50/60Hz (option 100, 110, 127V 50/60Hz)	230V-50/60Hz; 120V-50/60Hz (option 100, 110, 127V 50/60Hz)	230V-50/60Hz; 120V-50/60Hz (option 100, 110, 127V 50/60Hz)
POWER [W]	120	190	230
CAPACITY [ML]	8 x 15; 2 x 50; 2 x 60 PRP	500, 4 x CYTO	100, (24 capillaries 75 mm)
MAX. RPM [M⁻¹]	300 ÷ 4000, step 100 RPM	90 ÷ 6 000, step 1 RPM	90 ÷ 18 000, step 1 RPM
MAX. RCF [X G]	2 469	4 830	24 270
SETTING TIME	1 ÷ 99 min, step 1 min, ∞ (HOLD)	1s ÷ 99h 59min 59s, ∞, step 1s	1s ÷ 99h 59min 59s, ∞, step 1s
DIMENSIONS H X W X D [MM]	300 x 365 x 450	299 x 357 x 451	299 x 357 x 451
HEIGHT WITH OPEN LID [MM]	572	572	572

"MPW MED. INSTRUMENTS" SPÓŁDZIELNIA PRACY

ul. Boremlowska 46, 04-347 Warszawa
 tel.: +48 22 610 56 21 / fax: +48 22 610 55 36
 mail: mpw@mpw.pl / www.mpw.pl

REGON: 000401029
 KRS: 0000128941
 VAT UE: PL 5250001753
 BDO 000042924

Alior Bank S.A.
 16 2490 0005 0000 4530 5256 8362 **PLN**
 PL 44 2490 0005 0000 4600 2441 7250 **EUR**

WEIGHT [KG]	~ 21/23*	~ 22/23,7*	~ 20/21,5*
MIN. ACCELERATION TIME [S]	35	5	14
MIN. BRAKING TIME [S]	33	11	12
	* weight of the centrifuge depends on the voltage version	* weight of the centrifuge depends on the voltage version	* weight of the centrifuge depends on the voltage version

"MPW MED. INSTRUMENTS" SPÓŁDZIELNIA PRACY

ul. Boremlowska 46, 04-347 Warszawa
 tel.: +48 22 610 56 21 / fax: +48 22 610 55 36
 mail: mpw@mpw.pl / www.mpw.pl

REGON: 000401029
 KRS: 0000128941
 VAT UE: PL 5250001753
 BDO 000042924

Alior Bank S.A.
 16 2490 0005 0000 4530 5256 8362 **PLN**
 PL 44 2490 0005 0000 4600 2441 7250 **EUR**