

M-UNIVERSAL

M-DIAGNOSTIC

M-SCIENCE



PDF



v2





Easier access to the centrifuge chamber made of stainless steel.

#### FEATURES

modern design

low noise level

maintenance-free induction motor

modern software system

large graphic LCD - blanking the screen after a period of inactivity

2 screen modes (simplified and normal)

10 menu languages (PL, EN, ES, IT, PT, DE, RU, FR, SE, CZ)

visual and sound alerts of different notifications

100 user programs

10 acceleration/deceleration curves

speed or RCF adjustment

manual centrifugal radius adjustment (with automatic RCF correction)

time calculation from start key pressing or from reaching preselected speed

ascending or descending time counting

continuous operation mode - HOLD

SHORT time operation mode

changing of parameters during centrifuging

automatic rotor recognition

automatic lid opening

selected functions blocking, password protection

effective ventilation system

stainless steel rotor chamber

#### SAFETY

unbalance sensor

lid locking during rotor running

start blocking at opened lid

emergency lid lock release

product conforming with the EN-61010-1 and EN-61010-2-020, EN-61010-2-101 standards

# WYROBY MEDYCZNE DLA PROFESJONALISTÓW MEDICAL DEVICES FOR PROFESSIONAL USE ONLY



PN-EN ISO 9001  
PN-EN ISO 13485



## PRODUCENT WIRÓWEK LABORATORYNYCH MANUFACTURER OF LABORATORY CENTRIFUGES

### MPW MED. INSTRUMENTS

ul. Boremlowska 46 | 46 Boremlowska Street  
04-347 Warszawa, Polska | 04-347 Warsaw, Polska

[www.mpw.pl](http://www.mpw.pl)