


# PORTFOLIO



EN



## S-LINE

## M-LINE



**MPW-54**

**MPW-55**

**MPW-56**

**MPW-150R**

**M-BASIC**

**M-UNIVERSAL**



120 ml

48 ml

120 ml

90 ml

120 ml / PRP

500 ml

RCF  
max.

1 137,3 122 x g

15 279 x g

3 341 x g

21 382 x g

2 469 x g

24 270 x g

RPM  
max.

3 500, 5 800

14 500

6 000

15 000

4 000

18 000



-

-

-

-20°C ÷ +40°C

-

-



We have a great variety of equipment for each type of centrifuge which provides users with extensive research and diagnostic capabilities.

Products offered by MPW allow efficient and failure-free operation for a long time.



## M-LINE

## L-LINE



### M-DIAGNOSTIC

### M-SCIENCE

### MPW-260R

### MPW-260RH

### MPW-351e

### MPW-352

500 ml

100 ml

500 ml

500 ml

800 ml

1 000 ml

4 830 x g

24 270 x g

24 270 x g

24 270 x g

3 600 x g

24 270 x g

6 000

18 000

18 000

18 000

4 500

18 000

-

-

-20°C ÷ +40°C

-20°C ÷ +55°C

-

-



We offer various solutions for specialized laboratories.

We take up challenges to meet the requirements of cytodiagnosics by offering cytology kit to detect cellular changes in body fluids.

**PLEASE NOTE.** We reserve the right to change information included in our promotional materials at any time in the case the online offer changes. MPW MED. INSTRUMENTS is not responsible for printing mistakes.

## L-LINE



**MPW-352R**

**MPW-352RH**

1 000 ml

1 000 ml

30 065 x g

30 065 x g

18 000

18 000

-20°C ÷ +40°C

-20°C ÷ +55°C



## XL-LINE



**MPW-380**

**MPW-380R**

3 000 ml

3 000 ml

31 876 x g

31 876 x g


18 000

18 000

-


-20°C ÷ +40°C





## CYTO SET

base



microscope side




filter card



insert



PP tube 2,2ml  
with a cap



# MEDICAL DEVICES FOR PROFESSIONAL USE ONLY



PN-EN ISO 9001  
PN-EN ISO 13485



## MANUFACTURER OF LABORATORY CENTRIFUGES

### MPW MED. INSTRUMENTS

ul. Boremlowska 46 | 46 Boremlowska Street  
04-347 Warszawa, Polska | 04-347 Warsaw, Polska

[www.mpw.pl](http://www.mpw.pl)