


PORTFOLIO




EN



S-LINE

M-LINE



MPW-54

MPW-55

MPW-56

MPW-150R

M-BASIC

M-UNIVERSAL



120 ml

48 ml

120 ml

90 ml

120 ml / PRP

500 ml

RCF
max.

1 137,3 122 x g

15 279 x g

3 341 x g

21 382 x g

2 469 x g

24 270 x g

RPM
max.

3 500, 5 800

14 500

6 000

15 000

4 000

18 000



-

-

-

-20°C ÷ +40°C

-

-



We have a great variety of equipment for each type of centrifuge which provides users with extensive research and diagnostic capabilities.

Products offered by MPW allow efficient and failure-free operation for a long time.



M-LINE

L-LINE



M-DIAGNOSTIC

M-SCIENCE

MPW-260R

MPW-260RH

MPW-351e

MPW-352

500 ml

100 ml

500 ml

500 ml

800 ml

1 000 ml

4 427 x g

24 270 x g

24 270 x g

24 270 x g

3 600 x g

24 270 x g

6 000

18 000

18 000

18 000

4 500

18 000

-

-

-20°C ÷ +40°C

-20°C ÷ +55°C

-

-



We offer various solutions for specialized laboratories.

We take up challenges to meet the requirements of cytodiagnosics by offering cytology kit to detect cellular changes in body fluids.

PLEASE NOTE. We reserve the right to change information included in our promotional materials at any time in the case the online offer changes. MPW MED. INSTRUMENTS is not responsible for printing mistakes.

L-LINE



MPW-352R

MPW-352RH

1 000 ml

1 000 ml

30 065 x g

30 065 x g

18 000

18 000

-20°C ÷ +40°C

-20°C ÷ +55°C



XL-LINE



MPW-380

MPW-380R

3 000 ml

3 000 ml

31 876 x g

31 876 x g

18 000

18 000

-


-20°C ÷ +40°C






CYTOSET

base



microscope side




filter card



insert



PP tube 2,2ml
with a cap



WYROBY MEDYCZNE DLA PROFESJONALISTÓW MEDICAL DEVICES FOR PROFESSIONAL USE ONLY



PN-EN ISO 9001
PN-EN ISO 13485



PRODUCENT WIRÓWEK LABORATORYNYCH MANUFACTURER OF LABORATORY CENTRIFUGES

MPW MED. INSTRUMENTS

ul. Boremlowska 46 | 46 Boremlowska Street
04-347 Warszawa, Polska | 04-347 Warsaw, Polska

www.mpw.pl