


NEW ROTOR 12900  
for MPW-380 and 380R  
4 X 750 ml



PDF



v2

rotor	RPM [max]	RCF [max]	R [max]	angle	bucket	adaptor	tube cat. No	tube type	amount of tubes per rotor
REF: 12900 RPM max.: 4 100 RCF max.: 3 458 x g R max.: 184 mm Angle: 90°	4100	3458	184	90	13929	-	*	750 ml bottle [100 x 130 mm], Herolab®	4
	4100	3458	184	90		14901	*	15 ml tube, conical bottom, with cap [17 x 120 mm], Falcon®	56
						15050		15 ml tube, conical bottom, with cap [17 x 120 mm], Sarstedt®	56
	4100	3458	184	90			*	50 ml tube, conical bottom, with cap [30 x 117 mm], Falcon®	28
						14902	15052	50 ml tube, conical bottom, with cap [30 x 117 mm], Sarstedt®	28
	4100	3458	184	90			*	50 ml tube, conical bottom, without skirt [30 x 115 mm], Greiner®	28
	4100	3458	184	90			15048	15 ml Thermo Nalgene® [16 x 113 mm]	88
	4100	3458	184	90			15118	10 ml glass tube [16 x 100 mm]	88
	4100	3458	184	90			*	BD Vacutainer® [16 x 100 mm], [2,5-11 ml]	88
	4100	3458	184	90			*	Greiner Vacuette® [16 x 100 mm], [7-9 ml]	88
	4100	3458	184	90			*	10 ml Thermo Nalgene® Oak Ridge [16 x 81,5 mm]	88
	4100	3458	184	90			15119	7 ml glass tube [12 x 100 mm]	100
	4100	3458	184	90			15120	5 ml glass tube [12 x 75 mm]	100
	4100	3458	184	90			*	BD Vacutainer® [13 x 75 mm], [1,6-5,3 ml]	100
	4100	3458	184	90			*	BD Vacutainer® [13 x 100 mm], [4-7 ml]	100
	4100	3458	184	90			14904	Greiner Vacuette® [13 x 75 mm], [1-4,5 ml]	100
	4100	3458	184	90			*	Greiner Vacuette® [13 x 100 mm], [3,5-6 ml]	100
	4100	3458	184	90			*	Sarstedt S-Monovette® [13 x 65 mm], [2,6; 2,9; 3,4; 3,8 ml]	100
	4100	3458	184	90			*	Sarstedt S-Monovette® [13 x 75 mm], [2,7; 3; 4,3 ml]	100
	4100	3458	184	90			*	Sarstedt S-Monovette® [13 x 90 mm], [4,9; 5,6 ml]	100

\* tube is outside MPW offer



ul. Boremlowska 46  
04-347 Warszawa Polska  
46 Boremlowska Street  
04-347 Warsaw Poland

**MPW MED. INSTRUMENTS**  
www.mpw.pl

E-mail: mpw@mpw.pl  
Phone: +48 22 610 56 67  
Phone: +48 22 740 20 25  
Fax: +48 22 610 81 01

