

M-BASIC



HYBRID SWING-OUT AND ANGLE ROTOR

90°
/ 45°

BASIC



PROGRAM RPM/RCF $\times 100$

TIME

COVER

PROGRAM

SAVE

1
2
3
4
5

2 40

12 0



RCF

SHORT

COVER



FEATURES

horizontal swing-out and angle rotor [unique system for placing buckets in horizontal and angular positions]

maintenance-free induction motor

low noise level

5 user programs [chosen with use of direct selection key]

changing of parameters during centrifuging

saving previously used setting

clear LED display

device's status indicated by visual and sound signals

easy speed or RCF adjustment [chosen with use of direct selection key]

automatic lid opening

time calculation from START key pressing or from reaching preselected speed

2 acceleration characteristics: FAST, SOFT

4 deceleration characteristics: FAST, SOFT, LONG, RUN

HOLD operation mode [infinite time of centrifugation]

SHORT time operation mode

centrifugation cycles count

effective ventilation system

SAFETY


stainless steel rotor chamber

imbalance detection

lid locking during rotor running

start blocking at opened lid

product conforming with the EN-61010-1
and EN-61010-2-020, EN-61010-2-101 standards



WYROBY MEDYCZNE DLA PROFESJONALISTÓW MEDICAL DEVICES FOR PROFESSIONAL USE ONLY



PN-EN ISO 9001
PN-EN ISO 13485



PRODUCENT WIRÓWEK LABORATORYNYCH
MANUFACTURER OF LABORATORY CENTRIFUGES

MPW MED. INSTRUMENTS

ul. Boremlowska 46 | 46 Boremlowska Street
04-347 Warszawa, Polska | 04-347 Warsaw, Polska

www.mpw.pl